

Required Refrigerator and Freezer Specifications

November 2024

Overview

All NCIP participants must ensure their vaccine storage units meet the following specifications to guarantee the proper maintenance of vaccine temperature conditions. These units must be stand-alone refrigerators and freezers dedicated solely to vaccine storage.

Storage Unit Requirements

Key Requirements for Vaccine Storage Units:

- **Dedicated Use:** Storage units must be used exclusively for vaccines.
- **Temperature Range Compliance:**
 - **Refrigerators:** Must maintain a constant temperature between 2°C and 8°C (36°F to 46°F) for refrigerated vaccines.
 - **Freezers:** Must maintain a constant temperature between -15°C and -50°C (5°F to -58°F) for frozen vaccines.
- **Sufficient Capacity:** Units must have enough space to store the largest anticipated vaccine inventory for the year, without overcrowding. This includes space for water bottles in the refrigerator and frozen coolant packs in the freezer to stabilize temperatures.
- **Year-Round Temperature Reliability:** Units must be capable of reliably maintaining the required temperatures throughout the year.
- **Location:** Storage units should be placed in a well-ventilated room with at least 4 inches of space from the wall.
- **Warning Signs:** Post clear "DO NOT UNPLUG" signs at electrical outlets, on storage units, and circuit breakers to prevent accidental disconnection by staff, custodians, or maintenance personnel.

Refrigerator Requirements

The preferred refrigerator for vaccine storage is a stand-alone pharmaceutical grade or purpose-built unit designed for the precise storage of sensitive biologics. These refrigerators are the most secure option for vaccine storage.

Pharmacy-grade under the counter units and full-sized units are allowed for vaccine storage.

Key characteristics of the ideal unit include:

- **Consistent Temperature Maintenance:** Must reliably maintain temperatures between 2°C and 8°C (36°F to 46°F).
- **Fan Forced Air Circulation:** Ensures uniform temperature throughout the unit.
- **Microprocessor or Digital Thermostat:** Preferred for greater accuracy and stability.
- **Stable Performance:** Designed to handle ambient temperature fluctuations without compromising vaccine storage conditions.
- **Probe Port:** To allow monitoring equipment to be safely connected.

Freezer Requirements

The preferred freezer for vaccine storage is a stand-alone pharmaceutical-grade or purpose-built unit designed for the precise storage of sensitive biologics. These freezers are engineered to maintain stable temperatures required for vaccine preservation.

Types of Acceptable Freezers:

- **Manual Defrost Freezers:** These freezers require manual intervention to remove frost buildup but can effectively maintain the necessary temperature range.
- **Auto-Defrost (Frost-Free) Freezers:** These freezers automatically defrost at regular intervals. While this feature can be convenient, it can sometimes lead to temperature fluctuations during the defrost cycle. **Verify with the manufacturer** that the freezer will maintain the required temperature range (between -15°C and -50°C or 5°F to -58°F) even during the defrost cycle.
 - If purchasing a frost-free freezer, ensure that it can maintain in-range temperatures even when defrosting, as some models may temporarily exceed the -15°C minimum during these cycles.
 - A unit advertised as **-20°C capable** should be able to consistently maintain temperatures between **-15°C and -50°C**, including during defrost cycles.

Determining Freezer Size:

The **size of the freezer** required depends on the **enrollment population** and how frequently the NCIP participant orders vaccines. The following **Ordering Tiers** help determine the appropriate freezer size based on ordering frequency:

Table 1.1 Vaccine Ordering Tier and Frequency Guide

Ordering Tier	
Doses Per Year	Order Frequency
6,000 doses or more	Monthly
800 to 5,999 doses	Bi-monthly (every other month)
200 to 799 doses	Quarterly (every three months)
199 doses or less	As Needed

Purchasing New Vaccine Storage Units

Step 1: Assess Storage Needs

- The first step in purchasing a freezer is assessing the **clinical space available** and understanding the **workflow characteristics** of the practice.
- Larger Practices:** Should consider a **main storage area** with smaller units placed in the clinical area for easier access.

Step 2: Choose the Right Unit

- Stand-alone freezers can be purchased from **local retail stores** that sell general refrigeration units.
- Pharmacy-grade units** or purpose-built freezer models may take longer to acquire, but they are specifically designed for storing sensitive biologics like vaccines. These models often require service by **commercial refrigeration technicians**.

Step 3: Check Warranty & Service Options

- Be sure to **check the warranty** of the freezer before purchasing, as it may include restrictions on approved service providers. Most pharmacy-grade freezers should be serviced by **certified refrigeration technicians**.

Vaccine Storage Capacity and Equipment Guide

Table 2.1 Refrigerated and Freezer Vaccine Storage Capacity

Total Clients	Refrigerator Stand alone or pharmacy grade	Maximum Doses Refrigerator	Freezer Stand alone or pharmacy grade	Maximum Doses Freezer
Less than 50 Clients	4 to 6 cubic ft	Up to 200 doses	1.5 to 1.8 cubic ft	Up to 100 doses
Less than 300 Clients	8 to 11 cubic ft	200 to 400 doses	2 to 2.5 cubic ft	Up to 200 doses
Less than 700 Clients	15 to 25 cubic ft	400 to 900 doses	3 to 7 cubic ft	201 to 500 doses
700 + Clients	Multiple units	900 to 2,000 doses	7.1 to 15 cubic ft	501 to 6,000 doses

Table 2.2 Recommended Refrigerator Models with Electronic Microprocessors (2°C to 8°C)

Refrigerator product name/model numbers	Manufacturer Website
Accucold Pharmacy & Vaccine Refrigerators	www.accucold.com
AccuVax Vaccine Management System	trumedsystems.com/products/accuvax/
American Biotech Pharmacy Refrigerator Series	americanbiotechsupply.com
CliniCool Silver Series Pharmacy/Vaccine Refrigerators	www.labrepro.com

Migali Scientific Pharmacy/Vaccine Storage Refrigerator Evolution & Genesis Series	www.migaliscientific.com
Powers Scientific, Inc. Lab Single-Door, Pharmacy & 2-8°C Constant Temperature Refrigerators	powersscientific.com
So-Low Laboratory & Pharmacy 2°C to 8°C Refrigerator Models & Chromatography Refrigerator Models	so-low.com
TemPure Scientific Refrigerator 2°C to 8°C Models	tempurescientific.com
Thermo Scientific TSX Series High-Performance Pharmacy Refrigerators	thermoscientific.com

Table 2.3 Freezer models that maintain an optimum of at least -20°C:

Freezer product name/model numbers	Manufacturer/Supplier Website
Accucold Pharmacy & Vaccine Freezers	www.accucold.com
AccuVax Vaccine Management System	trumedsystems.com/products/accuvax/
CliniCool Silver Series Pharmacy/Vaccine Freezers	labrepco.com
EdgeStar Freezer 1.1 Cubic ft. CMF151L-1	edgestar.com
Futura Silver Series Laboratory Freezers	labrepco.com
Helmer i.Series & Horizon Series Laboratory Freezers	helmerinc.com
So-Low Laboratory Freezers	so-low.com
Thermo Scientific Forma High-Performance Auto Defrost Lab Freezers, Jewett High-Performance Lab Freezers, PL6500 Lab Freezers, ES Series Lab Freezers	thermoscientific.com

Table 2.4 Additional Manufacturer and Supplier Information

Manufacturer/Supplier	Website	Manufacturer/Supplier	Website
Accucold	accucold.com	Marvel Scientific	marvelscientific.com
Aegis Scientific	aegisfridge.com	Migali Scientific	migaliscientific.com
American Biotech	americanbiotechsupply.com	PHCbi	phchd.com/apac/biomedical
Compact Appliance	compactappliance.com	So-Low	so-low.com
Fisher Scientific	fishersci.com	Summit Appliance	summitappliance.com
Helmer	helmerinc.com	Thermo Scientific	thermoscientific.com

K2 Scientific	k2sci.com	VWR International	vwr.com
---------------	-----------	-------------------	---------

Note: The manufacturers and suppliers listed in *Table 2.4* offer units that are preset to specific temperatures suitable for vaccine storage. NCIP does not recommend a particular supplier but encourages you to compare options for the best fit and value for your facility. Prior to purchase, ensure that the unit will reliably maintain the temperature range required for vaccine storage: **Refrigerators:** 2°C to 8°C (with a set point between 4°C or 5°C)/**Freezers:** -20°C (-4°F), with a maximum temperature of -15°C (5°F)

For assistance, contact the NCIR Help Desk by
 Phone: 1.877.USE.NCIR (873-6247)
 Email: ncirhelp@dhhs.nc.gov